

Applications

- ATTO science
- Time-resolved spectroscopy
- Frontend for high power ultrafast amplifier systems
- Femtochemistry
- Coherent THz generation
- Materials processing
- OPA pumping


FEMTOPOWER™

compact™ PRO CE-Phase



FEMTOPOWER™ compact™ PRO CE-Phase is an ultrafast Ti:Sapphire multipass amplifier using a sealed vacuum chamber system for the amplifier crystal and Dispersive Mirror (DM) technology for dispersion compensation. The system generates the shortest, phase stabilized, high energy optical pulses commercially available with unprecedented quality, stability, and reproducibility. This compact multipass system is the optimum light source for attosecond science.

The amplifier is equipped with a CE-Phase stabilization system. It includes the FEMTOSOURCE™ rainbow™ DFG as seed source. The feedback signal for CE-Phase stabilization is derived from a new compact and stable DFG (Difference Frequency Generation) scheme utilizing the > 300 nm bandwidth (-10 dB) of the FEMTOSOURCE™ rainbow™ (*Optics Letters*, vol. 30, no. 3, Feb. 2005).



CE-Phase data (inloop) from a FEMTOLASERS CEP-stabilized FEMTOPOWER™ compact™ PRO.

FEMTOPOWER™

compact™ PRO CE-Phase

Extraordinary Features

- Excellent short & long term CEP stability
- Highest pulse quality and stability
- Broad spectrum | Ultrashort pulses
- Excellent pulse contrast
- Ultrahigh peak power
- Ultrafast oscillator output

Ultrashort Pulses

By employing an ultra broadband multipass scheme and DM technology, high quality, sub-30 fs pulses (sub-25 fs pulses on request) can be generated from a Ti:Sapphire kHz amplifier constituting the state of the art of high power high repetition rate lasers. (*Optics Letters, vol. 22, No. 20, pp. 1562, Oct. 15, 1997*)

Efficiency & Pulse quality

The use of DM technology allows the elimination of diffraction gratings, giving the highest efficiency and best pulse quality available from any commercial high power ultrafast amplifier.

Stability & Reliability


The amplifier crystal is cooled in a vacuum chamber, reducing thermal lensing and preventing dust mediated damage to the amplifier crystal. The compact footprint, as well as the superior mechanical and optical design of the entire system yields long term stability. Sufficient saturation of the Ti:Sapphire guarantees high pulse-to-pulse stability.

Inherent minimum CE-Phase noise guaranteed by:




- Dispersive Mirrors
- Large bandwidth
- CE-Phase decoupled from beam position
- Monolithic setup
- PPLN based
- Enhanced S/N ratio
- Gratingless
- Min. dispersion
- Min. CE-Phase Distortion


FP compact™ PRO CEP - SIDE VIEW, Dimensions in [mm] ([in])



FP compact™ PRO CEP - FRONT VIEW, Dimensions in [mm] ([in])



FP compact™ PRO CEP - TOP VIEW, Dimensions in [mm] ([in])



FEMTOPOWER™ compact™ PRO	Amplifier CEP	Amplifier CEP HP/HR	Oscillator (optional exit)
Pulse duration	< 30 fs < 25 fs	< 30 fs < 25 fs	< 7 fs CEP stabilised
Bandwidth (FWHM) @ 800 nm	> 40 nm > 60 nm	> 40 nm > 60 nm	> 250 nm @ -10 dB
Repetition rate	1 kHz	3 kHz 4 kHz 5 kHz	78 MHz
Output energy	> 800 μJ	> 800 μJ	> 2.1 nJ
Peak power	> 30 GW	> 30 GW	> 300 kW
Beam diameter (1/e ²) Beam divergence	15 mm (nom.) < 2 mrad	15 mm (nom.) < 2 mrad	< 2 mm < 2 mrad
Beam pointing stability (optional active stab.)	< 10 μrad rms @ Ø 15 mm	< 10 μrad rms @ Ø 15 mm	-
Spatial mode	TEM ₀₀ (M ² < 2)	TEM ₀₀ (M ² < 2)	-
Polarization	linear, horizontal	linear, horizontal	linear, horizontal
Contrast ratio ^{1) 2)}	> 10 ⁸ :1 or > 10 ⁴ :1	> 10 ⁸ :1 or > 10 ⁴ :1	-
Pulse to pulse energy stability	< 1.5 % rms	< 1.5 % rms	-

1) Defined as ratio between output pulse peak power and power of amplified spontaneous emission (ASE)
 2) Defined as ratio between output pulse peak power and power of any pulse > 1 ns before or after output pulse.

All specifications are subject to change without notice



VISIBLE AND INVISIBLE LASER RADIATION
 AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION

LASER TYPE: Ti:Sapphire
 WAVELENGTH: 650 - 1050 nm
 PULSE DURATION: > 15 fs
 MAX. PULSE ENERGY: 2 mJ @ 10 kHz
CLASS 4 LASER PRODUCT

compact™ PRO



FEMTOLASERS Produktions GmbH
 Fernkorngasse 10 | 1100 Wien | Austria
 P: +43 1 503 7002 0
 F: +43 1 503 7002 99
 info@femtolasers.com

FEMTOLASERS, Inc.
 124 Mt. Auburn St. | Suite 200 North
 Cambridge | MA 02138 | USA
 P: +1 978 456 9920
 F: +1 978 456 9922
 info@femtolasers.com



FEMTOLASERS' laser products are certified to comply with the Federal Regulations (21 CFR Subchapter J) as administered by Center of Devices and Radiological Health on all systems ordered for shipment after October 1, 2003.