



FEMTOSOURCE™
rainbow™

The **NEW** standard
for ultrafast oscillators

May 12th, 2010

The first CEP stabilized ultrafast oscillator with multiple parallel output options centered at 800 nm and 1030 nm. **Ultrafast** sub-7 fs pulse duration with 200 mW average power. **Ultrabroad** 260 nm spectral bandwidth with 200 mW for seeding applications. **Ultraclean** 2.5 nJ for spectra up to 50 nm in bandwidth for high contrast seeding applications. **Ultrapower** nearly 1 MW of peak power (500 mW, 7 fs pulses, non-CEP). **Ultra`red`** 1030 nm, ideal for seeding OPCPA pump lasers.

FEMTOLASERS Produktions GmbH

is the premiere manufacturer of ultrafast, compact and reliable laser oscillator and amplifier solutions. Founded in 1994 in Vienna, Austria, the company offers products and excellent services which have evolved to be the first choice among OEMs and scientists worldwide for demanding applications and solutions. FEMTOLASERS state-of-the-art products offer technology of highest quality.

We LIVE in an ULTRAFST world™

join us !